MIPARTM

--Spotlight 安装激活教程--

一、安装

注: MIPAR 软件要求版本: V4.5.0 版本以上(如有需要,请联系科普乐工程师)。

- 1. 安装包下载地址: Spotlight 安装包下载
- 2. "双击"打开安装包,默认 C 盘。点击 "Next"下一步。

😕 Setup - MIPAR Spotlight	-		\times
Select Destination Location Where should MIPAR Spotlight be installed?			
Setup will install MIPAR Spotlight into the following folder. To continue, click Next. If you would like to select a different folder, click Browse.			
C:\Program Files\MIPAR	Br	owse	
At least 7.57 GB of free disk space is required.			
Next	>	Ca	ncel

3. 点击"Install"安装。

😢 Setup - MIPAR Spotlight — 🗌	×	😕 Setup - MIPAR Spotlight — 🗆 🗙
Ready to Install Setup is now ready to begin installing MIPAR Spotlight on your computer. Click Install to continue with the installation, or click Back if you want to review or change any settings.		Installing Please wait while Setup installs MIPAR Spotlight on your computer.
Destination location: C:\Program FilesMIPAR		C \Program Files\MIPAR\plugina\Spotlight\spotlight_accurate.pt
< Back Install	Cancel	Cancel

4. 点击"Finish"完成。



二、激活

注: 要求 MIPAR 软件 V4.5.0 及以上版本; 安装好 Spotlight;下载提供的 V2C 许可证 (如有需要,请联系科普乐工程师)。

1. 双击安装的 MIPAR 软件。

2. 打开任意模块,点击"Help" -> "Activate License" -> "Update"。

3. 单击 "Activate",选择邮箱下载的 "Spotlight 的 V2C" 文件。

4. 重新启动 MIPAR 软件, Spotlight 和 Snap 已经安装激活完毕。加载图形即可 正常使用。



注: 一台电脑不能试用第二次, 如有问题, 请联系客服: 024-31014540